

stryker

Stryker's
titanium implants
inspired by biology



Why titanium?

Titanium nails have long been preferred for intramedullary nail fixation.* However, the use of titanium for other orthopaedic implants, especially in long bone plating, has become a much discussed topic in recent years. Titanium implants behave more like bone than stainless steel implants, offering a variety of advantages discussed in the following pages.

While orthopaedic implants are generally composed of either titanium or stainless steel, there are many differences in types of titanium - from grade 1 commercially pure titanium (CP-Ti) to grade 5 titanium alloy. There are also differences in the types of grade 5 titanium alloy, with Stryker utilizing the alloy Ti6Al4V for many years, including in the original Gamma and T2 nailing systems.



*Based on internal sales data.

Support fracture healing


Lower modulus of elasticity in titanium plates may support bone healing

Construct stiffness has been shown to directly correlate with secondary fracture healing.^{1,2}

Titanium's modulus of elasticity is closer to that of bone than stainless steel.³ Additionally, the design of the plate combined with the flexibility of titanium may allow the plate to more easily conform to the bone, allowing the surgeon to achieve a closer plate-bone fit.

Commercially pure titanium (CP-Ti) grade 2 and Ti6Al4V have a modulus of elasticity of 100-110 GPa and 100-130 GPa respectively, while 316L stainless steel has a Young's modulus of 200 GPa.⁴ The greater elasticity of titanium allows for more implant flexibility which directly affects callus formation.^{1,3} Interfragmentary motion in the millimeter range is proven to induce bone healing that does not occur with less motion.^{5,6} In comminuted fractures treated with bridge plating, the fracture must heal with external callus formation.¹

In a study by Rodriguez et al, there were 36 nonunions (13.3%) in a series of 271



Rodriguez EK, Boulton C, Weaver MJ, Herder LM, Morgan JH, Chacko AT, Appleton PT, Zurakowski D, Vrahas MS. Predictive factors of distal femoral fracture nonunion after lateral locked plating: a retrospective multicenter case-control study of 283 fractures. *Injury* Mar 2014, Vol 45, 554-559.

supracondylar fractures across 3 level 1 urban centers.^{2,*} While the authors were unable to eliminate design as a factor, it is worth noting that of the 239 fractures treated with titanium plates, 9.6% resulted in nonunions.

Of the 32 fractures treated with stainless steel, 40.6% resulted in nonunions.^{**}

The authors were able to interpret that both plate design and material, combined, had an impact.

*Study not specific to Stryker plating systems.

**"Nonunion was defined as the need for a secondary surgical procedure to improve poor healing (ie, bone grafting, implant exchange, other)."

*** Based on a predictive algorithm for nonunion using maximum likelihood estimation in logistic regression. Limitations of the study included that it was retrospective and there was no consistency around defining nonunions.

Answer

important clinical questions



Titanium implants may provide a clearer CT image, allowing you to more accurately answer important clinical questions

Titanium implants provide a clearer CT image than other common metals used for orthopaedic hardware, resulting in decreased image artifacts.^{8,9} Artifacts caused by metallic hardware may inhibit the ability to answer clinical questions.⁸



Metal objects in the field of view during a CT scan may reduce the effect of x-rays. "When an image is reconstructed using this incomplete data, it leads to unnatural changes in appearance, known as artifacts, which are often observed as

bright or dark streaks in the image."¹⁰

As shown in the pictures below, a study showed less artifact in titanium than stainless steel due to stainless steel's higher fluoroscopic density.⁸



Effect of metal composition on metal-related artifacts at multidetector CT. **(a)** Scout view of an ex vivo model of a pig femur with a 4.5-mm titanium alloy screw (arrow) and a 4.5-mm stainless steel screw (arrowhead). **(b)** CT image obtained through the short axis of the screws, with 140 kVp, 120 mAs, 0.6-mm collimation, and 1.0-mm reconstructed section thickness, shows more severe streak artifacts from the stainless steel screw than from the titanium alloy screw.¹⁸



Effect of peak voltage on metal-related artifacts. A comparison of CT images obtained through the short axis of the screws in the pig femur model, with 140 kVp **(c)** and 80 kVp **(d)** and with other parameters the same (120 mAs, 0.6-mm collimation, and 1.0-mm reconstructed section thickness), shows a less obtrusive artifact on the image obtained with a higher peak voltage.¹⁸

Increased tissue biocompatibility

Titanium implants are well known for their excellent biocompatibility, resulting in a favorable tissue response³

Titanium is well known for its excellent biocompatibility. When implanted into the body, titanium and its alloys remain essentially unchanged due to their excellent corrosion resistance.¹¹

For those patients with a nickel allergy, estimated to affect 17% of women and 3% of men, titanium implants are preferred over stainless steel.¹²

Although not shown to correlate with clinical data, internal in vitro studies show a 20.7% reduction in bacterial colonization on anodized Ti6Al4V samples compared to stainless steel.¹³

Titanium biocompatibility is defined by:¹¹



Titanium plating

reengineered

Ease of removal in long bone plating

Titanium plating has come a long way from its initial introduction. Competitive, early generation titanium locking plates were impacted by screw stripping and cold-welding, making plate removal very difficult.¹⁴

Since the 2013 launch of AxSOS 3 Titanium, only 0.001% of cases have reported concerns with screw removal most often due to the overpowering of screws.¹⁵


During the development of AxSOS 3, the plates were implanted into sheep tibiae and extracted at 6 months and 12

months. This preclinical study demonstrated that all screws were extractable when following the surgical technique and manually tightening locking screws using a torque limiter.¹⁶

Design features built into the AxSOS 3 screw/plate interface

to allow for consistent screw insertion and removal¹⁶

- 1. Patented, round threaded screw head design** – instead of the typical sharp threaded design used in the shaft of screws, the rounded threads on the locking screw head make it easier to tighten and release the locking screw from the plate¹⁶
- 2. Conically shaped screw head** – makes it easier for the screw to find its path to help ensure good axial alignment¹⁶
- 3. One-threaded screw in a dual-threaded plate** – one thread path is clear¹⁶
- 4. Type II anodization** – a surface treatment that partially melts the surface of the implant, changing the chemical composition, and resulting in up to a 17% increase in fatigue strength and a reduction in friction and wear^{16,17}



References:

1. Henderson CE, et al. 2010 Mid-America Orthopedic Association Physician in Training Award: Healing Complications are Common After Locked Plating for Distal Femurs. *Clinical Orthopaedics and Related Research* (2011) 469:1757-1765.
2. Rodriguez EK et al. Mechanical construct characteristics predisposing to non-union after locked lateral plating of distal femur fractures. *Journal of Orthopaedic Trauma* 30(8), August 2016.
3. Pohler O EM. Unalloyed titanium for implants in bone surgery. *Injury* 31(2000) S-D7-13.
4. Kulkarni M et al. Nanomedicine, Chapter: Biomaterial Surface Modification Of Titanium and Titanium Alloys for Medical Applications, Publisher: One Central Press, Editors: Alexander Seifalian, Achala de Mel, Deepak M. Kalaskar, pp.111-136
5. Claes LE, et al. Effects of mechanical factors on the fracture healing process. *Clin Orthop Relat Res.* 1998;355(suppl): S132-S147.
6. Kenwright J, et al. Effect of controlled axial micro-movement on healing of tibial fractures. *Lancet.* 1986;2:1185-1187.
7. Rodriguez EK, Boulton C, Weaver MJ, Herder LM, Morgan JH, Chacko AT, Appleton PT, Zurakowski D, Vrahas MS. Predictive factors of distal femoral fracture nonunion after lateral locked plating: a retrospective multicenter case-control study of 283 fractures. *Injury* Mar 2014, Volume 45, pp. 554-559.
8. Kataoka ML et al. A review of factors that affect artifact from metallic hardware on multi-row detector computed tomography. *Curr Probl Diagn Radiol*, July/August 2010 (pp. 125-136).
9. Berg BV et al. Multi-detector CT imaging in the postoperative orthopedic patient with metal hardware. *Eur J Radiol* 2006;60:470-9.
10. Gjestebj L et al. Metal artifact reduction in CT: where are we after four decades? *IEEE Access*, pp. 5826-5949, October 2016.
11. Ratner D. "A Perspective on Titanium Biocompatibility." *Titanium in Medicine: Material Science, Surface Science, Engineering, Biological Responses and Medical Applications.* Eds. Brunette D M, et al. Heidelberg:Springer-Verlag, 2001, pp. 1-10.
12. Thyssen JP et al. Metal allergy – a review on exposures, penetration, genetics, prevalence, and clinical implications. *Chem Res Toxicol.* 2010 Feb 15;23(2):309-18.
13. Data on file. Test report STO050615MA.01F
14. Suzuki T et al. Technical problems and complications in the removal of the Less Invasive Stabilization System. *Journal of Orthopaedic Trauma*, 24(6), June 2010.
15. Data on file. As of March 18, 2016, USA only.
16. White paper, content ID# AXSOS-WP-4. Konowalczyk, S. et al. AxSOS 3 Titanium Plating System – Screw Removal. September 2014.
17. Data on file. Test report 300709CB1
18. Lee MJ, et al. Overcoming artifacts from metallic orthopedic implants at high-field-strength MR imaging and multi-detector CT. *RadioGraphics*, May-June 2007, Vol. 27, Number 3 pp. 791-803.

Plating that fits.

Fit for patients

Using the powerful functionality of SOMA, our population-based design environment, we are able to create anatomically shaped plates that are enhanced for a broad range of patients.

Fit for procedures

Intelligent instruments, smart design features and thoughtful innovation allow for a streamlined procedure, letting you focus on what really matters.

Fit for you

Dedicated and hard-working, our sales force provides you with market-leading support.

Fit for the future

Evolving our plating portfolio is a continuous process, and our ambition is to lead in all aspects of trauma plating.



Trauma & Extremities

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate a Stryker product. A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area. Stryker Corporation or its affiliates own, use, or have applied for the following trademarks or service marks: Anchorage, Asnis, AxSOS 3, Fixos, Gamma, Gamma3, Plating that fits., SOMA, Stryker, T2, VariAx. All other trademarks are trademarks of their respective owners or holders.

Manufacturer:

Stryker GmbH
Bohnackerweg 1
2545 Selzach
Switzerland
stryker.com



CE 0123