

Linear Reduction Clamp

Exposure. **Reduction.** Fixation.

Modernized Reduction.

Ratchet-based mechanism of action

- ▶ Allows for consistent and repeatable compressive forces

Four hook options

- ▶ Four different hook options with varying sizes and tip ends to suit your indication needs. Hook attachments can be oriented in **8W** different ways relative to the handle and can be changed intra-operatively.

Removable handle technology

- ▶ The handle can be removed from the assembly without any loss of compressive forces. This allows for multiple hook attachments to be assembled from the same set.


Instrumentation made simple.

Cannulated compression rod for direct screw

- ▶ Cannulated compression rod can accept $\varnothing 2.7\text{mm}$ cortical, $\varnothing 3.5\text{mm}$ cortical, $\varnothing 4.5\text{mm}$ cortical, and $\varnothing 6.5\text{mm}$ cancellous screws

Compatible with PRO system

- ▶ Complements our existing PRO pelvis portfolio



A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery. The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker. All other trademarks are trademarks of their respective owners or holders.

Content ID: XXXX, 09-2021
Copyright © 2021 Stryker



Manufactured by:
Avalign Technologies
8727 Clinton Park Drive
Fort Wayne, IN 46825
1.877.289.1096
www.avalign.com

Distributed by:
Stryker
stryker.com